



NO *Touch*™



THIS

NOT THIS

**CONTAMINATION AND TEARS
OF SURGICAL INSTRUMENT PACKS
ARE A THING OF THE PAST**


 **Material Handling USA**
800-326-4403
www.mh-usa.com




NO *Touch*TM

CONTAMINATION AND TEARS ARE A THING OF THE PAST WITH NO-TOUCHTM... FROM THE SPD TO THE OR, WITH VIRTUALLY NO HUMAN INTERACTION

The No Touch system eliminates human interaction during instrument transport thanks to our dedicated ISO Standard Stainless Steel shelves and baskets and our functional lines of open and closed transportations carts.




STEP 1
Clean Instrument
Wrapping



STEP 2
Wrapped Instruments are
placed on No Touch Shelf and
into Sterilizer (shelf is handled,
pack is not touched)



STEP 3
No Touch Shelf (w/sterile pack) is
transported to Sterile Storage (shelf is
handled, pack is not touched)



STEP 4
No Touch Shelf (w/sterile
pack) is transported to OR to
be Opened by Nurse


Mobile & Stationary Storage Systems



Aluminum Transport Cart



Open Transportation Trolley



Sterile Instrument Pack Storage

Eliminate contamination and tears during transportation by touching only the carrier shelf, and not the sterile pack

- ✓ Maximum sterility during transportation
- ✓ Safeguards patients from infection



NOT THIS



THIS


STAINLESS STEEL MODULES

Modules suitable for sterilization


Material Handling USA




60X40cm (24"X16") | 10cm (4") Deep
WS-B104060-SS



60X40cm (24"X16") | 20cm (8") Deep
WS-B204060-SS



Flat 60X40cm (24"X16")
SSS-6040WH



Flat 30X40cm (12"X16")
SSS-3040WH